

OCUNOVIS BIOHANCE Gel Eye Drops is formulated to provide lubrication and aqueous barrier where eyes may be unprotected or where the natural blink reflex is diminished. BIOHANCE is a unique, cross-linked hyaluronic acid (HA) that has been specifically modified for the ocular surface environment in dogs and cats. It creates a cellular matrix of modified HA that provides a sheer aqueous and lubricates the ocular surface.

Physical barrier product for topical application to animals.

FOR USE IN VETERINARY HOSPITALS.

For topical application to the eye where long lasting lubrication may be of benefit in dogs and cats.

FOR ANIMAL USE ONLY.

**Administration:**

1. Carefully remove cap from bottle.
2. Place 1 or 2 drops directly into the eye. The natural movement of the eyelid will distribute the drops over the eye surface.
3. OCULENIS BIOHANCE can be applied as little as 2 times daily or as eye lubrication is needed to per veterinarian recommendation.
4. To avoid contamination, do not touch the tip of the container to any surface.
5. Replace cap securely after using.

**Ingredients:**

All ingredients are naturally occurring in the body and do not cause pain or stringing when applied:

- Water
- BIOHANCE hyaluronic acid (0.4%)
- BioSalts:
  - Sodium Chloride
  - Disodium Phosphate
  - Potassium Chloride
  - Potassium Phosphate

Contains no preservatives or excipients that may irritate the eyes. Free of antibiotics or steroids

**Disposal:**

Dispose of container by wrapping with paper and putting in the garbage.

**Storage:**

Store below 30°C (room temperature). May be stored in the refrigerator. If refrigerated, the product should be allowed to warm to room temperature prior to use. DO NOT FREEZE. Store upright. KEEP OUT OF REACH OF CHILDREN.

Exclusively distributed by:  
Provet PTY LTD  
133a Lenore Drive  
Erskine Park NSW 2756

Manufactured in USA by SentrX Animal Care, Inc.

Net Contents: 5mL